

SOP 99-01 Generating immortalized macaque B-lymphoblastoid cell lines

Purpose and principle:

Autologous cells are used at targets in assays to detect cytotoxic T cell activity. This SOP describes the generation of immortalized B-LCL from macaque monkey PBMC using a transforming Herpesvirus.

Materials:

- Macaque whole blood
- RPMI 1640 with 12% FCS
- *Herpesvirus papio* containing supernatant (from S594 cell cultures, for monkey B-LCL)
- 15 ml polypropylene tubes
- 15 ml polystyrene tubes
- Pipets
- 96-well plates
- Table-top centrifuge
- 25 cm² flasks
- 75 cm² flasks

Procedure:

1. Ficoll isolate lymphocytes from monkey blood and lyse any red cells.
2. Plate 10⁵ lymphocytes per well in a 96 well plate in 100 µl of culture medium.
3. Add 100 µl of S594 supernatant or 100 µl media for controls.
4. Incubate at 37°C, 5% CO₂, changing 100 µl of the medium at least once a week.
5. Dark foci of immortalized cells should appear within 2-4 weeks.
6. Expand when possible and maintain cell lines at 0.4x10⁶ cells/ml in 12% FCS RPMI in T25 flasks.

Special handling/waste disposal:

Treat *H.papio* as a biohazard and handle at BSL2+.

Aruna Seth
Beth Israel Deaconess Medical Center
Date: 19 September 1999
Revised: